

MIXER MILL MM 500 CONTROL



The MM 500 control is a high energy laboratory ball mill that can be used for dry, wet and cryogenic grinding with a frequency of up to 30 Hz. It is the first mixer mill in the market that allows to monitor and control the temperature of a grinding process.

The temperature area covers a range from -100 to 100 °C. For maximum flexibility, the mill can be operated with various thermal fluids, enabling the use of different tempering devices for cooling or heating. If liquid nitrogen is chosen for cooling, the mill needs to be equipped with the optionally available extension device cryoPad. The innovative cryoPad technology allows to select and control a specific cooling temperature in the range from -100 to 0 °C for the grinding process.



[Click to view video](#)

Product Video

MIXER MILL MM 500 VARIO

Variability redefined



PROCESSING 2 TO 50 SAMPLES IN ONE BATCH

The new Mixer Mill MM 500 vario is a versatile bench-top unit which provides ultimate performance with maximum flexibility for your sample preparation process.

It is used for dry, wet and cryogenic grinding of small sample amounts with high throughput. The MM 500 vario can be equipped with screw-top grinding jars from 1.5 ml to 50 ml. Available materials include hardened steel, stainless steel, tungsten carbide, agate, zirconium oxide, PTFE.

For biological applications such as homogenization of plant materials, tissues or for cell disruption via bead beating, the MM 500 vario can be equipped with different adapters for single-use vials from 0.2 - 5 ml.



[Click to view video](#)

Product Video

PERFORMANCE AND DESIGN

- | Grind sizes down to 5 µm possible
- | Powerful pulverization with up to 35 Hz
- | Twice as fast as mills with only 30 Hz
- | Ergonomic design with touch display

FLEXIBILITY

- | Suitable for dry, wet or cryogenic grinding
- | Equally suited for rapid (<2 min) and long-term grinding (up to 99 hours)
- | 12 SOPs and 4 program cycles with max. 99 repeats facilitate routine applications
- | Grinding jars of various sizes and materials; adapters for single-use vials



SOLUTIONS FOR BIOLOGICAL APPLICATIONS AND CELL DISRUPTION

Mixer mills are widely used for homogenizing biological samples such as tissue, liver, muscle, plants, corn or sputum. For cell disruption via bead beating mixer mills are also the perfect solution. The MM 500 vario accommodates adapters for different single-use vials:

0.2 ml / 1.5 ml / 2 ml / 5 ml

Depending on the tube and the type of adapter, 5-10 tubes can be inserted per adapter, resulting in the following maximum vial capacity per batch:

50 x 1.5 or 2 ml (or 60 x 0.2 ml)

20 x 5 ml


Typical grinding / homogenizing processes of biological materials take less than 2 min. Cell disruption with a good reproducibility and efficiency takes 20 sec to 5 min, depending on the cell type. Cryogenic grinding in single-use vials can also be a suitable way to pulverize tough or temperature-sensitive samples.



Solutions for biological applications and cell disruption

MAXIMIZE YOUR THROUGHPUT

- | Maximum frequency of 35 Hz reduces grinding time to 50% compared to mixer mills with only 30 Hz
- | Six grinding stations instead of two, as in classic mixer mills, ensure high throughput
- | In combination, the throughput of the MM 500 vario is 6 times higher compared to a classic mixer mill



Grinding of basalt in the MM 500 vario results in better fineness compared to classic Mixer Mills thanks to the increased frequency of 35 Hz instead of 30 Hz (50 ml jar + 12 x 12 mm grinding balls).

FOR SAFE AND EFFECTIVE GRINDING PROCESSES
ACCESSORIES FOR THE MM 500 VARIO



GRINDING JARS IN 6 DIFFERENT MATERIALS

All grinding jars of the smaller MM 400 model are compatible with the MM 500 vario. The nominal volume ranges from 1.5 ml to 50 ml and available materials include hardened steel, stainless steel, agate, tungsten carbide, zirconium oxide and PTFE, ensuring contamination-free sample preparation.




ADAPTERS FOR SINGLE-USE VIALS

The adapters for 0.5 / 1.5 / 2 / 5 ml single-use vials of the smaller MM 400 model can be used in the MM 500 vario.

TYPICAL SAMPLE MATERIALS

RETSCH mixer mills are true allrounders. They homogenize, for example, waste, soil, chemical products, coated tablets, drugs, ores, grain, tissue, glass, hair, ceramics, bones, plastics, alloys, minerals, oil seeds, plants, sewage sludge, pills, textiles, wool etc.



concrete



wood



plastic granulates



hair

To find the best solution for your sample preparation task, visit our application database.

TECHNICAL DATA

Applications	mechanochemistry, homogenization, size reduction, cryogenic grinding, mechanical alloying, mixing,
Field of application	agriculture, biology, chemistry / plastics, construction materials, engineering / electronics, environment / recycling, food, geology / metallurgy, glass / ceramics, medicine / pharmaceuticals
Feed material	hard, medium-hard, soft, brittle, elastic, fibrous
Size reduction principle	impact, friction
Material feed size*	<= 8 mm
Final fineness*	~ 5 µm
Batch size / feed quantity*	max. 6 x 20 ml
No. of grinding stations	6
Setting of vibrational frequency	digital, 3 - 35 Hz (180 - 2100 min ⁻¹)
Typical mean grinding time	30 s - 2 min
Dry grinding	yes
Wet grinding	yes
Cryogenic grinding	yes
Cell disruption with reaction vials	yes
Self-centering clamping device	yes
Type of grinding jars	screw top design
Material of grinding tools	hardened steel, stainless steel, tungsten carbide, zirconium oxide, agate, PTFE
Grinding jar sizes	1.5 ml / 5 ml / 10 ml / 25 ml / 35 ml / 50 ml
Setting of grinding time	digital, 10 s - 8 h
Total grinding time	99 h
Storable SOPs	12
Number of storable cycle programs	4 (with 99 repeats)
Electrical supply data	100-120V, 50/60 Hz; 200-230V, 50/60Hz
Power connection	1-phase
Protection code	IP 30
Power consumption	750 W
W x H x D closed	690 x 375 x 585 mm
Net weight	~ 60 kg
Standards	CE

*depending on feed material and instrument configuration/settings

MIXER MILL MM 500 VARIO

FUNCTIONAL PRINCIPLE

The grinding jars of the mixer mill MM 500 vario perform radial oscillations in a horizontal position. The inertia of the grinding balls causes them to impact with high energy on the sample material at the rounded ends of the grinding jars and pulverize it. Also, the movement of the grinding jars combined with the movement of the balls result in the intensive mixing of the sample. The degree of mixing can be increased even further by using several smaller balls.



[Click to view video](#)

THE NEW RETSCH APP




As the leading solution provider for sample preparation equipment, RETSCH has taken operating convenience to the next level and created the new RETSCH App. This tool makes working with your RETSCH mill easy and convenient:

- | Operate your devices via your smart phone or tablet
- | Control your devices based on your own application routines
- | Access information from the RETSCH database
- | Get in touch with the RETSCH service team via the app

www.retsch.com/mm500-vario

MIXER MILL MM 500 NANO

Comfortable production of particles in the nanometer range



MORE THAN AN ALTERNATIVE TO A PLANETARY BALL MILL

The mixer mill MM 500 nano is a compact, versatile bench-top unit which has been developed specially for dry, wet and cryogenic grinding of up to 2 x 45 ml sample material within seconds. With a maximum frequency of 35 Hz, it generates enough energy to produce particles in the nanometer range.

The robust high-performance drive makes the mill suitable for long-term grinding processes up to 99 hours and thus very interesting for research and mechanochemistry.

Hence, the MM 500 nano is a unique mixer mill in the market to provide a real alternative to grinding in a planetary ball mill – with more comfortable handling and less warming effects.



[Click to view video](#)

ADVANTAGES THROUGH DESIGN

- | Very easy, comfortable clamping and handling of the grinding jars
- | Jars can stay clamped while taking a sub-sample or visual checks of fineness
- | Ergonomic design with touch display for easy parameter setting
- | 12 SOPs & 4 program cycles with up to 99 repeats to facilitate routine applications

FLEXIBILITY

- | Equally suited for rapid pulverization <2 min and long-term grinding up to 99 hours
- | Use one large grinding ball in the High Impact mode or several smaller balls in the High Friction mode
- | Use the MM 500 nano for routine sample preparation applications, for nano-grinding or for research applications such as mechanochemistry and mechanical alloying


MIXER MILL MM 500 NANO

GRIND SIZES IN THE NANOMETER RANGE

- | New jar design allows for optimized usage of jar volume also for wet grinding
- | Final fineness < 100 nm possible thanks to maximized energy input at 35 Hz
- | Less warming effects, thus grinding can usually be done without grinding breaks for cooling down

Result: You get your nano sample within the shortest time.

NANO-GRINDING OF TITANIUM DIOXIDE WITHOUT COOLING BREAKS



Nano grinding of 25 g titanium dioxide in a 125 ml grinding jar zirconium oxide with 275 g balls 0.1 mm, 30 ml 1% NaPO₄ solution. A particle size of 90 nm was achieved after 120 minutes of grinding.

NEW SCREW-LOCK GRINDING JARS: BENEFITS FOR YOU

The new Screw-Lock grinding jars of the MM 500 nano are suitable for dry, wet or cryogenic grinding. The jar sizes (50 / 80 / 125 ml) are larger than those of the classic mixer mills, allowing for grinding 2 x 45 ml per batch. The jars are pressure-tight up to 5 bar, the integrated safety closure allows for convenient handling. The new jar design is very beneficial for wet grinding and pulverizing fibrous samples like hair.




Thanks to the flat lid, the nominal volume can be fully used, for instance when milling fibrous samples, or to ensure the optimum mixture of sample, small grinding balls and liquid for wet grinding.

Available materials include hardened steel, stainless steel, tungsten carbide and zirconium oxide ensuring contamination-free sample preparation. Aeration lids for all jar sizes and materials are available, e.g. for grinding under inert atmosphere.

MIXER MILL MM 500 NANO

TYPICAL SAMPLE MATERIALS

RETSCH mixer mills are true allrounders. They homogenize, for example: alloys, animal feed, bones, ceramics, chemical products, coal, coke, drugs, electronic scrap, glass, grains, hair, minerals, oil seeds, ores, paper, plant materials, plastics, sewage sludge, soils, straw, tablets, textiles, tissue, tobacco, waste samples, wood, wool, etc.



titanium oxide
wet grinding



metal alloy
dry grinding



hair
dry grinding




tyre rubber
cryogenic grinding

To find the best solution for your sample preparation task, visit our application database.

MIXER MILL MM 500 VARIO

NANO GRINDING OF ALUMINUM OXIDE IN THE MM 500 NANO

Narrow particle size distribution of aluminium oxide after grinding (Nano grinding of 30 g aluminium oxide in a 125 ml grinding jar zirconium oxide with 275 g balls 0.1 mm, 33 ml 0.5% NaPO₄ solution)



TECHNICAL DATA

Applications	mechanochemistry, mechanical alloying, size reduction, mixing, homogenization, cryogenic grinding
Field of application	agriculture, biology, chemistry / plastics, construction materials, engineering / electronics, environment / recycling, food, geology / metallurgy, glass / ceramics, medicine / pharmaceuticals
Feed material	hard, medium-hard, soft, brittle, elastic, fibrous
Size reduction principle	impact, friction
Material feed size*	<= 10 mm
Final fineness*	~ 0.1 µm
Batch size / feed quantity*	max. 2 x 45 ml
No. of grinding stations	2
Setting of vibrational frequency	digital, 3 - 35 Hz (180 - 2100 min ⁻¹)
Typical mean grinding time	30 s - 2 min
Dry grinding	yes
Wet grinding	yes
Cryogenic grinding	yes
Cell disruption with reaction vials	no
Type of grinding jars	screw-lock with integrated safety closure devices
Material of grinding tools	hardened steel, stainless steel, tungsten carbide, zirconium oxide
Grinding jar sizes	50 ml / 80 ml / 125 ml
Setting of grinding time	digital, 10 s - 8 h
Total grinding time	99 h
Storable SOPs	12
Number of storable cycle programs	4 (with 99 repeats)
Electrical supply data	100-120V, 50/60 Hz; 200-230V, 50/60Hz
Power connection	1-phase
Protection code	IP 30
Power consumption	750 W
W x H x D closed	690 x 375 x 585 mm
Net weight	~ 60 kg
Standards	CE

*depending on feed material and instrument configuration/settings

MIXER MILL MM 500 NANO

FUNCTIONAL PRINCIPLE

The grinding jars of the mixer mill MM 500 nano perform radial oscillations in a horizontal position. The inertia of the grinding balls causes them to impact with high energy on the sample material at the rounded ends of the grinding jars and pulverize it. Also, the movement of the grinding jars combined with the movement of the balls result in the intensive mixing of the sample. The degree of mixing can be increased even further by using several smaller balls.



[Click to view video](#)

THE NEW RETSCH APP



As the leading solution provider for sample preparation equipment, RETSCH has taken operating convenience to the next level and created the new RETSCH App. This tool makes working with your RETSCH mill easy and convenient:

- | Operate your devices via your smart phone or tablet
- | Control your devices based on your own application routines
- | Access information from the RETSCH database
- | Get in touch with the RETSCH service team via the app

www.retsch.com/mm500-nano


ADVANTAGES THROUGH DESIGN

- | Dry, wet and cryogenic grinding with up to 30 Hz for high energy grinding
- | Fast and comfortable sample processing with two screw lock jars of up to 125 ml each
- | Patented hermetically closed fluid system ensures the safe operation of thermal fluids
- | Wide range of accessories available, including ventilation lids and heavy-metal-free grinding jars (also for cryogenic grinding)
- | Ergonomic jar clamping, low noise level, user friendly parameter setting via touch display




TEMPERATURE MONITORING AND CONTROL

- | Continuous temperature monitoring throughout a grinding process
- | Cooling and heating in a range from -100 to 100°C
- | Operation is possible with liquid nitrogen or other thermal fluid
- | High flexibility in terms of selecting a tempering device for temperature regulation (LN₂ supply, cryostat, chiller, ...).
- | Low temperature grinding is possible without LN₂



CRYOPAD


- | Extension device cryoPad is required for the operation with LN₂
- | The cryoPad regulates the flow of LN₂ through the thermal plate
- | The cryoPad technology allows to select and maintain a specific cooling temperature in the range between -100 and 0°C while using LN₂



MIXER MILL MM 500 CONTROL

TEMPERATURE REGULATION BASED ON THERMAL PLATES

The cooling and heating of the sample material is realized with the patented concept of thermal plates, making sample cooling with, e. g., open liquid nitrogen baths or dry ice obsolete. For tempering, the grinding jars are simply placed on top of the thermal plates. When the grinding jars come in contact with the thermal plates, heat is effectively transferred from or to the jars via the tempering device. The patented hermetically sealed fluid design allows to operate the mill with different thermal fluids, ensuring a flexible and safe temperature regulation and requiring only minimal effort for the user. Depending on the operational setup that is built up, the temperature of the thermal plates can be set in the range from -100 to +100 °C.



MIXER MILL MM 500 CONTROL CONFIGURATIONS

To control the temperature of a grinding process, the mill needs to be connected to an external tempering device. Basically, there are two options:

1. Temperature regulation with liquid nitrogen

The mill is operated with liquid nitrogen and connected to a nitrogen tank. In this setup the mill must be extended with the optionally available extension device cryoPad. The patented PID (proportional–integral–derivative) system of the cryoPad controls the flow of liquid nitrogen and herewith the temperature of the thermal plates. In this setup, it is possible to select and maintain the temperature of the thermal plates at a specific value. The desired temperature is adjusted via the touch display and can be selected within a range from -100 to 0 °C, in steps of 10.


Setup 1: Extension device cryoPad and LN₂ tank for the operation with liquid nitrogen.



2. Cooling or heating with a liquid thermal fluid

In this setup, the mill can either be connected to a cryostat, to a chiller or to the water tap. The external tempering device regulates the corresponding thermal fluid to a defined temperature and the thermal fluid transfers this temperature to the thermal plates. As during a grinding process, a significant amount of heat may also develop inside the jar, the temperature of the thermal plates may be manipulated. To sum up, the actual temperature of the thermal plates depends on both, on the temperature of the thermal fluid and on grinding parameters, like frequency, time, jar volume, size of grinding balls. For a maximum control of the grinding process, the actual temperature of the thermal plates is continuously monitored in the touch display.

Setup 2: Operation with an external tempering device; e.g. water tap, chiller or thermostat.



MIXER MILL MM 500 CONTROL

APPLICATION EXAMPLES

The temperature regulation of the MM 500 control is especially designed for the processing of temperature-sensitive sample materials. Cooling or heating may have different objectives.

Cooling can be used for example:

- | Preserving temperature-sensitive analytes (like volatile substances or pharmaceutical and food ingredients)
- | Embrittlement
- | Wet grinding below room temperature
- | Mechanochemistry

Some applications are improved if the sample material is heated up during the process. Examples for heating are:

- | Paste making (in food industry)
- | Intensifying mechanochemical reactions

The required temperatures and the operational setup depend on the specific application.



[Click to view video](#)

PRESERVING TEMPERATURE-SENSITIVE ANALYTES

Some analytes are modified, destroyed or vaporized if the sample material gets too warm. If specific temperature levels are exceeded, the structure of, for example, proteins, pharmaceutical substances or food ingredients may be essentially changed.

By keeping the temperature at a moderate level throughout the grinding process, temperature-sensitive natural substances are physically preserved in their original state for analysis.



Grinding of coffee beans at low temperatures for natural substance analysis.

CRYOGENIC GRINDING

Temperatures below 0 °C are suitable for the embrittlement and homogenization of for example ductile or sticky food. If heavy-metal-free grinding is required, jars of zirconium oxide or tungsten carbide can be used.


If cooled down to -100 °C, it is also possible to successfully embrittle some polymers.



Fast milling of black Fluoro-carbon rubber (FKM) by embrittling the sample in two 125 ml jars at -100 °C.

WET GRINDING < 30 °C


If using a chiller, powerful wet grinding can be performed at 30 Hz and below room temperature without considering any cooling breaks.



Particle size and temperature development for TiO₂ in a wet grinding process with 30 Hz and 2 x 125 ml jars

MECHANOCHEMISTRY

By cooling the sample throughout a mechanochemical process, the formation of undesired derivatives can be prevented. It is also possible to apply some heating, for example to initiate chemical reactions and increase product yields.



Fast formation of ZIF-8 at 30 Hz. By keeping the temperature below 0 °C, the formation of non-porous Zeolitic Imidazolate Frameworks (ZIF-8) is inhibited.

MIXER MILL MM 500 CONTROL

ACCESSORIES FOR MAXIMUM FLEXIBILITY

The MM 500 control is equipped with two grinding stations for screw-lock jars, available in the materials stainless steel, hardened steel, zirconium oxide and wolfram carbide in different jar sizes up to 80 ml or 125 ml. Heavy-metal-free grinding is possible, also at low temperatures. Ventilation lids allow for ventilation of the jars and for working under inert atmospheres.




MIXER MILL MM 500 CONTROL

TYPICAL SAMPLE MATERIALS

As the MM 500 control can be used with or without cooling, the mill offers a wide variety of applications. It can be used to homogenize, for example, waste, soil, chemical products, coated tablets, drugs, ores, grain, tissue, glass, hair, ceramics, bones, plastics, alloys, minerals, oil seeds, plants, sewage sludge, pills, textiles, wool etc.



raisins



coated tablets



polystyrene



soil

To find the best solution for your sample preparation task, visit our application database.

TECHNICAL DATA

MIXER MILL MM 500 CONTROL

Applications	mechanochemistry, mechanical alloying, size reduction, mixing, homogenization, cryogenic grinding
Field of application	agriculture, biology, chemistry / plastics, construction materials, engineering / electronics, environment / recycling, food, geology / metallurgy, glass / ceramics, medicine / pharmaceuticals
Feed material	hard, medium-hard, soft, brittle, elastic, fibrous
Size reduction principle	impact, friction
Material feed size*	<= 10 mm
Final fineness*	~ 0.1 µm
Batch size / feed quantity*	max. 2 x 45 ml
Grinding chamber volume	max. 2 x 125 ml
No. of grinding stations	2
Setting of vibrational frequency	digital, 3 - 30 Hz (180 -1800 min ⁻¹)
Setting of temperature setpoint	digital, 0 ... -100 °C (only with cryoPad)
Setting of sample cooling time	digital, 0 ... 60 min (only with cryoPad)
Setting of grinding time	digital, 10 s - 8 h
Total grinding time	99 h
Storable SOPs	12
Number of storable cycle programs	4 (with 99 repeats)
Typical mean grinding time	30 s - 2 min
Dry grinding	yes
Wet grinding	yes
Cryogenic grinding	yes
Type of grinding jars	screw-lock with integrated safety closure devices
Material of grinding tools	hardened steel, stainless steel, tungsten carbide, zirconium oxide
Grinding jar sizes	50 ml / 80 ml / 125 ml
Electrical supply data	100-120V, 50/60 Hz; 200-230V, 50/60Hz
Power connection	1-phase

Protection code	IP 30
Power consumption	750 W
W x H x D closed	690 x 375 x 585 mm
W x H x D closed with cryoPad	690 x 485 x 585 mm
Net weight	~ 63 kg
Standards	CE
Connection thread size device input	G 1/4" (inner thread)
Connection thread size tubing set	G 3/8" (outer thread)
Permissible operating pressure cooling device (provided by customer)	0 ... 5 bar
typical pressure range of continuous cooling unit e.g. cryostat	1 ... 2 bar
permissible pressure range of LN2 supply	1.2 ... 1.4 bar
Permissible fluids	water, water-glycole mixture, thermal oil, liquid nitrogen
Thermal applications	embrittling, cooling, heating, temperature control
temperature range of fluids	+100 °C ... -196 °C
temperature range of cooling plates	+100 °C ... -100 °C

*depending on feed material and instrument configuration/settings

TECHNICAL DATA

CRYOPAD

Applications	cryogenic grinding with liquid nitrogen
Interface	RS-232 (MM 500 control)
Communication connection	via included connection cable
Power supply	via external power supply
Electrical supply data (input external power supply)	100-230V, 50/60 Hz
External power supply classification	Medical grade isolation level
Electrical supply data (input cryoPad)	24 V, 1 A
Accessories	LN2 Autofill 150L, LN2 Autofill 50L
LED status light	yes
W x H x D	670 x 110 x 590 mm
Net weight	~ 26 kg
Standards	CE
Connection thread size device input	G 1/4" (inner thread)
Connection thread size of stainless steel tubing adapter	UNF 3/4"
Permissible pressure range of LN2 supply	1.2 ...1.4 bar
Permissible fluids	Liquid nitrogen
Emissions	Liquid nitrogen gas, condensation
Connection	via included tubing set
Exhaust outlet	via included Exhaust adapter and aluminum corrugated tube
temperature range of fluids	-196 °C
temperature control algorithm	PID temperature control
Setting of temperature setpoint	digital, 0 ... -100 °C
Setting of sample cooling time	digital, 0 ... 60 min

MIXER MILL MM 500 CONTROL

FUNCTIONAL PRINCIPLE

The grinding jars of the mixer mill MM 500 control perform radial oscillations in a horizontal position. The inertia of the grinding balls causes them to impact with high energy on the sample material at the rounded ends of the grinding jars and pulverize it. High energy milling is possible by operating at high frequencies up to 30 Hz. The movement of the grinding jars combined with the movement of the balls further causes grinding effects due to friction and additionally result in an effective mixing of the sample. The degree of mixing can be increased even further by using several smaller balls.



[Click to view video](#)

www.retsch.com/mm500-control